User's Manual

ELRA Electrode

EN	Electrode for electrosurgical device	1 page
DA	Elektrode til elektrokirurgisk indretning	5 page
DU	Elektrode voor elektrochirurgisch apparaat	10 page
FR	Électrode pour dispositif électrochirurgical	15 page
GE	Elektrode für elektrochirurgische Geräte	20 page
GR	Ηλεκτρόδιο για ηλεκτροχειρουργική συσκευή	25 page
IT	Elettrodo per dispositivo elettrochirurgico	31 page
NO	Elektrode for elektrokirurgisk enhet	36 page
PO	Elétrodo para dispositivo electrocirúrgico	41 page
SP	Electrodo para dispositivo electroquirúrgico	46 page
SW	Elektrod för elektrokirurgisk anordning	51 page
TU	Elektrocerrahi cihazı için elektrot	56 page

User's manual

1. Introduction

The electrode for this electrosurgical device is used by connecting it to the RF lesion generator from STARmed Co., Ltd. which was designed to coagulate tissues. Tissues can be coagulated with the supply of radiofrequency energy from the electric cautery, which is connected to the generator for monitoring the impedance and temperature of tissue. The ELRA Electrode includes both active and passive electrodes. Radiofrequency energy flows between the electrodes.

2. Guide for use

Read all of the following information and also the provided user's generator guide before use. ELRA Electrode is designed to be used in conjunction with RF lesion generator from STARmed Co., Ltd. If the guidelines are not properly followed, damages can be caused from electricity or heat and the equipment may malfunction.

3. Intended Use

Sterile single-use RF coagulation electrode is intended for coagulation of tissue during Endoscopic, percutaneous, laparoscopic, and intraoperative surgical procedures.

4. Intended patient population

age	Newborn to geriatric
weight	>2.5Kg
gender	Male and Female
health	Use for patients with liver cancer, thyroid cancer, lung cancer, kidney cancer or other lesions.
Nationality	Multiple
Patient state	The patient is not the device user. Not relevant unless the patient is agitated.

Contraindication group	Patients with implantable pacemaker and implantable cardioverter/defibrillators.
------------------------	--

WARRANTY

STARmed Co., Ltd. warranties sufficient care for the design and manufacture of this device. We do not warrant any content that is not described in this manual. The handling, storage, and operation of the ELRA Electrode beyond what is described in the IFU may directly affect the results obtained from its use. STARmed's device under this warranty is limited to repair and replacement, and STARmed Co., Ltd. is not responsible for any economic expenditures caused by unexpected issues, significant loss or damage. STARmed Co., Ltd. is not transferring responsibility and duty related to this device to other parties, and will not take any responsibility regarding reuse, and usage of re-sterilized or expired products.

5. Warning

- Before using this product, check whether there is any damage to the packaging.
- The expiration date is three (3) years from the sterilization date. Use of the products after the expiration date is prohibited.
- This device is disposable and reuse is prohibited.
- This device is designed to be used by a qualified and trained Medical doctor.
- This device should be stored in a cool place without sunlight and humidity.
- Do not touch the electrode's active tip while the power supply is on.
- Do not touch the shaft of electrode with fingers or other tools (Insulator and components are made of metallic).
- Do not touch the electrode's hand piece/cable with metallic materials. Touching the cable with metallic materials may cause the patient or operator to get shocked, cause a fire, or have an accident with

injury.

- Do not use on patients with cardiac pacemakers or other active implants.
- Do not use near conductive materials.
- Like all other electrosurgical devices, do not use in the presence of flammable anesthetics, nitrous oxide or any other or flammable substances. Doing this may cause a fire in the electrosurgical device.
- Flammable gases or other materials with a high risk of flammability should be prohibited from being placed nearby the RF generator. If a flammable detergent has been used for cleaning or disinfection, the liquid should have evaporated before the start of the procedure.
- Cloth and gauze soaked in oxygen can possibly be ignited from sparks which occur from normal use of the device. Please use with caution.
- Do not activate the electrode while it is touching any metallic substances or tools. Doing this may cause an unexpected injury to the patient or operator.
- Do not touch metals or tools with the electrode while the power supply is on. Doing this may cause an unexpected injury to the patient or operator and/or damage to the electrode or other devices.
- When using the electrode, do not hold the tip of the electrode by forceps, do not pass the electrode through a hole of other devices, and do not scratch the surface on the electrode's tip by sharp objects. These actions might destroy the insulation of the electrode and/or damage the electrode. It might cause of unexpected injury to patient or user.
- Do not access to lesion without proper imaging guidance. It might be cause an unexpected injury to patients and/or damage to electrode.

6. Potential complications

- Burning by overheating of surgical device.
- Dangers from an inexperienced operator's use.


- Side effects or cross infection from reuse.
- Weakness of liver functions.
- Delayed bleeding in the operated body parts.
- Recurrence of cancer.
- Perforation.
- Respiratory depression or arrest.
- Symptoms after RFA treatment including: abdominal pain, fever, nausea, headache, right shoulder joint pain and chest discomfort

The Potential complications is not limited to the list above. Other complications may occur.

7. Before Use

- (1) Before using this product, check whether there is any damage to the packaging.
- (2) Be careful when taking out the sterilized electrode from the packaging to avoid contamination.

 **Caution:** Before inserting the electrode into the lesion, check if the exposed length of the electrode and other specifications are identical with the specifications on the product's labels.

8. During use

- (1) Connect the electrode to the RF lesion generator. Confirm the connection status.

 **Caution:** Before inserting the electrode into the lesion, check if the exposed length of the electrode and other specifications are identical with the specifications on the labels.

- (2) Check the location of the EUS or scope when

targeting the lesion for accuracy.

 **Caution** : Refer to the user's manual for warnings and how to use the EUS or scope.

(3) Insert the guide wire into the instrument channel and position it exactly on the lesion so that the electrode can be easily placed on the lesion.

 **Caution**: The maximum outer diameter of a guide wire for use with this product is 0.035 inch.

 **Caution**: Please refer to the operating method and the cautions for use in the user's manual for the relevant guide wire product.

(4) Make sure to place the electrode's tip exactly on the lesion by connecting it to the guide wire that is already placed on the lesion through the instrument channel of the EUS or scope.

 **Caution**: After placing the electrode on the lesion, move the electrode during the RF procedure to avoid ablating an unwanted zone.

(5) Perform the RF procedure using the RF lesion generator under proper operating conditions (output, temperature, and time).

(6) When the RF procedure is over, stop generating power output by pressing the STOP button on the RF lesion generator.

9. After use

(1) When the procedure is finished, turn off the POWER switch of the RF lesion generator.


(2) Separate the electrode from the instrument channel of the EUS or scope.

 **Caution**: Be careful that the patient and/or operator are not injured from the tip of the electrode when removing the electrode from the EUS or scope. Follow the user's manual when removing the EUS or scope from the patient.


(3) Dispose of the electrode according to the waste materials disposal procedure.

10. Storage

(1) Store at room temperature.

 **Caution**: This device is intended to be sold and used by the order of medical doctor according to the related medical device laws.

11. Action for product damage

 **Caution**: If any damage on product is visually noticed, do not use to prevent injury to patient and/or operator.

12. Expiration date : Three (3) years from sterilization

13. Symbol

Authorized Representative in the European Community

Consult instructions for use

Manufacturer

Catalogue number, Model name

Sterilized using ethylene oxide

Caution

Use-by date

Batch code

Do not re-use

Date of manufacture

Reuse Precaution Statement

The contents of this product have been sterilized with E.O. gas. If the sterilized package is damaged, do not use and contact STARmed Co., Ltd.

This product can be used only one time. Reuse, re-treatment or re-sterilization is not allowed.

If the electrode is reused, re-treated, or re-sterilized,

it may cause defects in the structural functions of the electrode and damage of the RF lesion generator which may lead to injury, disease, or death to the patient. Reuse, re-treatment, or re-sterilization can cause device contamination and infections among patients or from other diseases. Contamination of this device may be a cause of injury, disease, or death.

14. Descriptions

Common name: Electrode for electrosurgical device

Trade name: ELRA Electrode

Catalogue number, Model name: Refer to the product label

Lot number: Refer to the product label

Manufactured date: Corresponding year and date

Number of items: 1 set

Sterilization method: E.O. gas

Expiration date: Three (3) years from date of sterilization

Method of use: Refer to the manual

Intended use: Refer to the manual

Cautions in use: Refer to the manual

Method of storage: Store at room temperature

Manufacturer: STARmed Co., Ltd.

※ This product is a disposable medical device and reuse is not allowed.

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Üretici firma: STARmed Co, Ltd
STARmed Co, Ltd