Niti-S

Tracheobronchial Stent

Table of Contents

English.....	3
Français.....	7
Deutsch	12
Italiano	17
Español	21
Português	26
Türkçe.....	30
Русский	34
Čeština	39
Српски.....	43
Dansk.....	47
Ελληνικά.....	51
Lietuvių.....	56
Polski	60
Română	65
Svenska.....	69
Norsk	73
Български език	77
中國語.....	82

MRI Information

MR Conditional

Niti-S stent was determined to be MR-conditional.

Non-clinical testing demonstrated that the Niti-S stent is MR Conditional according to ASTM F2503. A patient with this device can be scanned safely, immediately after placement under the following Conditions:

Static Magnetic Field

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 720-Gauss/cm or less

MRI-Related Heating

In non-clinical testing, Niti-S stent produced the following temperature rises during MRI performed for 15-min of scanning (i.e., per pulse sequence) in 1.5-Tesla/64-MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numaris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner) and 3-Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems:

	<u>1.5-Tesla</u>	<u>3-Tesla</u>
MR system reported, whole body averaged SAR	2.9-W/kg	2.9-W/kg
Calorimetry measured values, whole body averaged SAR	2.1-W/kg	2.7-W/kg
Highest temperature change	+1.9°C	+2.4°C

These temperature changes will not pose a hazard to a human subject under the conditions indicated above.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position Niti-S stent. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 10-mm relative to the size and shape of this implant.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	7,695-mm ²	676-mm ²	7,984-mm ²	977-mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

User's Manual

1. Description

The Niti-S Tracheobronchial Stent consists of the implantable metallic stent and introducer system.

The stent is made of Nitinol wire. It is a flexible, fine mesh tubular prosthesis which has radiopaque markers on each end and at the center.

Model Name
Niti-S Tracheal Uncovered Stent
Niti-S Tracheal Covered Stent
Niti-S Bronchial Uncovered Stent
Niti-S Bronchial Covered Stent

Figure 1. Stent Model

The Stent is loaded in introducer system and upon deployment the stent imparts an outward radial force on the luminal surface of the trachea and/or bronchus to establish patency.

Niti-S Full Covered Tracheobronchial Stents used in benign stricture can be removed; (see Warnings).

Niti-S Full Covered Tracheobronchial stents can be repositioned after deployment (see Warnings).

Stent deployment shall be performed under OTW (Over The Wire) procedure using a bronchoscope. Introducer systems are divided into two types of 12Fr or less and 14Fr or more (See Figure 2, 3).

Please note that the TTS (Through The Scope) approach shall not be used.


Figure 2. Introducer System (12Fr. or less)


Figure 3. Introducer System (14Fr. or more)

2. Principle of Operation

The outer sheath is pulled back by immobilizing the hub in one hand, grasping the Y-connector with the other hand, and gently sliding the Y-connector along the pusher (2nd inner catheter) towards the hub. Retraction of the outer sheath releases the stent.

3. Indication for Use

The Niti-S Tracheobronchial Stent is intended for maintaining Tracheal or bronchial luminal patency in tracheobronchial stricture caused by intrinsic and/or extrinsic malignant and/or benign stricture.

WARRANTY

Taewoong Medical Co., LTD. warrants that reasonable care has been applied within the design and subsequent manufacturing process of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Taewoong's control directly affect the instrument and the results obtained from its use. Taewoong's obligation under this warranty is limited to the repair or replacement of this instrument and Taewoong shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. Taewoong neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. Taewoong assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.

4. Contraindication

The Niti-S Tracheobronchial Stent is contraindicated for, but is not limited to;

- Placement in polypoid lesions.
- Patient with bleeding disorder.
- Intra-abdominal abscess/perforation.

- Patients with coagulopathy
- Strictures that do not allow passage of a guidewire.
- Any use other than those specifically outlined under indications for use.
- Removal or repositioning of fully deployed uncovered Stents is contraindicated. (See Warnings).
- Suspected or impending perforation.
- Recapturing a stent during its deployment is contraindicated.

5. Warnings

- The device should be used with caution and only after careful consideration in patients with elevated bleeding times, coagulopathies, or in patients with radiation colitis or proctitis.
- Chemoradiation therapy or radiotherapy alone may lead to tumor shrinkage and subsequent stent migration.
- The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.
- Do not expose the introducer system to organic solvent (e.g. Alcohol)
- Do not use with Ethiodol or Lipiodol contrast media.
- Niti-S Full Covered Stents can be removed where the treating doctor's clinical assessment of the stented stricture indicates the Stent can be safely removed. Caution should be exercised in deciding to and when removing the Stent.
- Niti-S Full Covered Stent cannot be removed when there is tumor in-growth/overgrowth/occlusion of the Stent lumen.
- Niti-S Full Covered Stent can be repositioned immediately after deployment.
- Uncovered/bare Stents should not be removed once fully deployed, (See Contraindications).
- Do not attempt to recapture/reload a stent once its deployment is advanced.

6. Potential complications

Potential complications associated with the use of Niti-S Stent may include, but are not limited to:

Procedural Complications

- Bleeding
- Stent misplacement
- Pain
- Death (other than that due to normal disease progression)
- Perforation
- Deployment failure

Post Stent Placement Complications

- Bleeding
- Pain
- Perforation
- Stent migration
- Stent occlusion
- Tumor overgrowth
- Tumor ingrowth
- Death (other than that due to normal disease progression)
- Fever
- Foreign body sensation
- Sepsis
- Infection
- Edema
- Pneumothorax
- Halitosis
- Hemoptysis
- Restenosis
- Dyspnea
- Stent fracture
- Hypoxia
- Cough
- Granulation tissue formation
- Stent misshaping
- Tracheoesophageal fistula
- Pneumonia
- Respiratory failure
- Retentions of secretions
- Atelectasis
- Ischemia

7. Equipment required

- Fluoroscope and/or Bronchoscope
- 0.035" / 0.89 mm guidewire (12Fr. or less)
- 0.038" / 0.97 mm guidewire (14Fr. or more)
- Introducer sheath appropriately sized for stent and introducer system

8. Precautions

Read the entire User's Manual thoroughly before using this device. It should only be used by or under the supervision of physicians thoroughly trained in the placement of stents. A thorough understanding of the techniques, principles, clinical applications and risks associated with this procedure is necessary before using the device.

- Care should be taken when removing the introducer system and guidewire immediately after stent deployment since this may result in stent dislodgement if the stent has not been adequately deployed.
- Care should be taken when performing dilation after the Stent has been deployed as this may result in perforation, bleeding, Stent dislodgement or Stent migration.
- The packaging and the device should be inspected prior to use.
- Use of fluoroscopy is recommended to ensure correct placement of the device.
- Check the expiration date "Use by". Do not use the device beyond the use by date.
- The Niti-S Stent is supplied sterile. Do not use any the packaging is opened or damaged.
- The Niti-S Stent is intended for single use only. Do not resterilize and/or reuse the device.

9. Instructions in the event of Damage

WARNING: Visually inspect the system for any sign of damage. DO NOT USE if the system has any visible signs of damage. Failure to observe this precaution may result in patient injury.

10. Procedure

① Examine stricture bronchoscopically and fluoroscopically

- a) Carefully examine both the proximal and distal segment of stricture bronchoscopically and/or fluoroscopically.
- b) The Internal luminal diameter should be measured exactly with bronchoscope and/or fluoroscope.

② Stent Size Determination

- a) Measure the length of the target stricture.
- b) Select a stent size that is 20 to 40mm longer than the measured length of the stricture in order to cover fully both ends of the lesion.
- c) Measure the diameter of the reference stricture - it is necessary to select a stent which has an unconstrained diameter about 1 to 4mm larger than the largest reference target diameter, to achieve secure placement.

③ Stent Deployment Preparation

- The Niti-S Stent can be placed with the aid of fluoroscopy, and/or bronchoscopy.
- Pass a 0.035" (0.89 mm) guidewire to the level of the stricture.
- Pass a 0.038" (0.97 mm) guidewire to the level of the stricture.

A. 12Fr. or less introducer system

- a) Under the bronchoscopic and/or fluoroscopic guidance, insert a bronchoscope to the level of the obstruction, then introduce the guide wire through the working channel of the bronchoscope. Advance until the guide wire across the target stricture to where the stent introducer system will be placed over the guide wire.
- b) Remove the stylet from the distal end of the introducer.
- c) Ensure that the valve of Y-connector connecting the inner sheath and outer sheath is locked by rotation proximal valve end in a clockwise direction to prevent premature stent deployment.
- d) Flush the inner lumen of introducer system.

B. 14Fr. or more introducer system

- a) Remove the stylet from the distal end of the introducer.
- b) Under bronchoscopy and/or fluoroscopy guidance, insert a guide wire across the stricture to where the stent introducer system will be placed over the guide wire.
- c) Ensure that the valve of Y-connector connecting the inner sheath and outer sheath is locked by rotation proximal valve end in a clockwise direction to prevent premature stent deployment.

④ Stent Deployment Procedure


Figure 4

PRECAUTION: Do not twist introducer system or employ a boring motion during the deployment as this may affect positioning and ultimate function of stent

- a) Remove the bronchoscope
- b) Under the fluoroscope and/or bronchoscope guidance, position the introducer system to the center of the target stricture exactly.
- c) Once the introducer system is in the correct position for deployment, unlock the proximal valve of the Y-connector by turning the valve more than twice in an anti-clockwise direction.

- d) To begin stent deployment, immobilize the hub in one hand and grasp the Y-connector with the other hand. Gently slide the Y-connector back along the pusher towards the hub.
- e) When the center X-ray marker reaches the center of target stricture, continue pulling back on the Y-connector until the stent is fully deployed. (See figure 4, 5)


Figure 5

CAUTION Do not push forward or pull backward on the hub with the stent partially deployed. The hub must be securely immobilized. Inadvertent movement of the hub may cause misalignment of the stent and possible damage to the target or stricture.

5 After Stent Deployment

- a) Examine the stent fluoroscopically and/or bronchoscopically to confirm expansion.
- b) Carefully remove the introducer system, guidewire and bronchoscope from the patient. If excessive resistance is felt during removal, wait 3~5 minutes to allow further stent expansion. (Place the inner sheath back into the outer sheath as the original state prior to removal.)
- c) Balloon dilatation inside the stent can be performed on demand.

11. Perform routine post implant procedure

- a) Assess the size and stricture of the Stent lumen. A Stent may require up to 1 to 3 days to expand fully.
- b) Doctor's experience and discretion can determine the appropriate drug regimen for each patient.
- c) After implantation, patient should remain on a soft diet until otherwise determined by the treating doctor.
- d) Observe the patient for development of any complications.

12. Instructions for removal of Niti-S Full Covered Stents (see Warnings)

Visually examine the Stent for any tumor ingrowth/overgrowth into the Stent lumen or whether the Stent is occluded. If the Stent lumen is clear, carefully remove using a forcep and/or snare. Grasp the retrieval string and/or collapse the proximal end of the Stent then carefully retrieve the Stent. If the Stent cannot be easily withdrawn, do not remove the Stent.

Caution: Do not allow excessive force to remove the stent as it may cause disconnect to the retrieval string.

To reposition a Niti-S Full Covered Stent immediately after deployment, use forceps or a snare to grasp the retrieval string and gently adjust to the correct placement.

Please note: the stent can only be repositioned and/or removed proximally.

Reuse Precaution Statement

Contents supplied STERILE (ethylene oxide (EO)). Do not use if sterile barrier is damaged. In the event of damaged packaging, call your Taewoong Medical Co., Ltd. representative. For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and /or cause patient infection or cross infection, including, but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Storage: Store at room temperature (10~40°C).

Disposal Requirements: The introducer system of Niti-S Tracheobronchial Stent must be properly sealed and disposed in compliance with the regulation of local or hospital at the end of its use.

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