UVENTA™

Ureteral Stent

Table of Contents

English.....	3
Français.....	7
Deutsch	11
Italiano.....	15
Español	19
Português	23
Türkçe.....	27
Čeština.....	31
Dansk.....	35
Ελληνικά.....	39
Lietuvių.....	43
Polski.....	47
Română.....	51
Svenska.....	55
Norsk.....	59

MRI Information

MR Conditional

UVENTA™ Ureteral Stent was determined to be MR-conditional.

Non-clinical testing demonstrated that the UVENTA™ Ureteral Stent is MR Conditional according to ASTM F2503. A patient with this device can be scanned safely, immediately after placement under the following Conditions:

Static Magnetic Field

- Static magnetic field of 3-Tesla or less

- Maximum spatial gradient magnetic field of 720-Gauss/cm or less

MRI-Related Heating

In non-clinical testing, UVENTA™ Ureteral Stent produced the following temperature rises during MRI performed for 15-min of scanning (i.e., per pulse sequence) in 1.5-Tesla/64-MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numaris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner) and 3-Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems:

	<u>1.5-Tesla</u>	<u>3-Tesla</u>
MR system reported, whole body averaged SAR	2.9-W/kg	2.9-W/kg
Calorimetry measured values, whole body averaged SAR	2.1-W/kg	2.7-W/kg
Highest temperature change	+1.9°C	+2.4°C

These temperature changes will not pose a hazard to a human subject under the conditions indicated above.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position UVENTA™ Ureteral Stent. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 10-mm relative to the size and shape of this implant.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	7,695-mm ²	676-mm ²	7,984-mm ²	977-mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

1. Description

UVENTA™ Ureteral Stent consists of an implantable and self-expandable metallic Stent and its delivery system.

UVENTA™ Ureteral Stent has a unique dual metal mesh structure formed by two stents and a membrane: Outer bare metal stent, inner bare metal stent, and PTFE membrane in between the two stents. The stent is composed of a flexible nitinol wire that forms a fine mesh tubular prosthesis and X-ray markers on it. The radiopaque markers help targeting the position of deployed stent.

Model Name
UVENTA™ Ureteral Stent

The Stent is pre-loaded in its delivery system and upon deployment the Stent imparts an outward radial force on the luminal surface of the Ureter to establish patency.


Figure1. Introducer system

2. Principle of Operation

The outer sheath is pulled back by immobilizing the hub in one hand, grasping the Y-connector with the other hand, and gently sliding the Y-connector along the pusher (2nd inner catheter) towards the hub. Retraction of the outer sheath releases the stent.

3. Indications for Use

UVENTA™ Ureteral Stent is intended for maintaining Ureteral luminal patency in Ureteral stricture caused by intrinsic and/or extrinsic malignant stricture.

WARRANTY

Taewoong Medical Co., LTD. warrants that reasonable care has been applied within the design and subsequent manufacturing process of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Taewoong's control directly affect the instrument and the results obtained from its use. Taewoong's obligation under this warranty is limited to the repair or replacement of this instrument and Taewoong shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. Taewoong neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. Taewoong assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.

4. Contraindication

The contraindications of UVENTA™ Ureteral Stent are as followed, but are not limited to:

- Bladder outlet obstruction
- Coagulation defect
- Untreated infection
- Spastic or noncompliant bladder
- Bladder fistula
- Recapturing a stent during its deployment is contraindicated.

5. Warnings

- The device should be used with caution and only after careful consideration in patients with elevated bleeding times, coagulopathies, or in patients with radiation colitis or proctitis.
- Chemoradiation therapy or radiotherapy alone may lead to tumor shrinkage and subsequent Stent migration.
- The Stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.
- Do not expose the delivery system to organic solvent (e.g. Alcohol)
- Do not use with Ethiodol or Lipiodol contrast media.
- UVENTA™ Ureteral Stent can be removed only when the charge surgeon's clinical assessment of the Stented stricture indicates the Stent can be safely removed. Cautions should be well-trained in decision of how and when to remove the Stent.
- Do not attempt to recapture/reload a stent once its deployment is advanced.

6. Potential complications

Potential complications associated with the use of UVENTA™ Ureteral Stent may include, but are not limited to:

Procedural Complications

- Bleeding
- Stent misplacement or inadequate expansion
- Pain
- Death (other than that due to normal disease progression)
- Luminal Perforation

Post Stent Placement Complications

- Bleeding
- Transient pain
- Perforation
- Urine reflux
- Urine incontinence
- Stent migration
- Stent re-occlusion due to encrustation or calcification.
- Stent occlusion due to tumor over-growth / hyperplasia
- Fever
- Foreign body sensation
- Death (other than that due to normal disease progression)
- Infection
- Edema
- Fistula

7. Required Equipments

- C-Arm (Fluoroscopic vision) and Cystoscope and / or Ureteroscope
- 0.035"/ 0.89 mm guidewire
- Ureteral catheter
- Balloon dilator (Ø 6mm)

8. Precautions

- Read the entire User's Manual thoroughly before using this device. It should only be used by or under the supervision of surgeon thoroughly trained for the placement of Stents. A thorough understanding of the techniques, principles, clinical applications and risk associated with this procedure is required before using the device.
- Care should be taken when removing the delivery system and guidewire immediately after Stent deployment since this may result in Stent dislodgement if the Stent has not been adequately deployed.
- Care should be taken when performing dilation after the Stent has been deployed as this may result in perforation, bleeding, Stent dislodgement or Stent migration.
- The packaging should be inspected prior to use as well as the device.
- Use of C-Arm (fluoroscopy) is recommended to ensure a correct placement of the device.
- Check the expiration date "Use by". Do not use the device beyond the labeled use by date.
- UVENTA™ Ureteral Stent is supplied sterile. Do not use the device if the packaging is opened or damaged.
- UVENTA™ Ureteral Stent is intended for a single use only. Do not resterilize and/or reuse the device.

9. Instructions in the event of Damage

WARNING: Visually inspect the system for any sign of damage prior to use. DO NOT USE if the system has any visible signs of damage. Failure to observe this precaution may result in patient injury.

10. Procedure

NOTICE: UVENTA™ Ureteral Stent can be performed by urologist and interventional radiologist. Following procedures indicate only fundamental guidance of UVENTA™ Ureteral Stent.

UVENTA™ Ureteral Stent must be operated in a sterile room.

① Examine stricture endoscopically and fluoroscopically

Carefully examine both the proximal and distal segment of stricture endoscopically and/or fluoroscopically.

② Stent Size Determination

- a) Measure the length of the target stricture.
- b) Select a Stent size longer than the measured length of the stricture in order to cover both ends of the stricture securely.

③ Stent Deployment Preparation

- UVENTA™ Ureteral Stent can be placed with the aid of fluoroscopy, and cysto / ureteroscopy.
- Advance the 0.035" (0.89 mm) guidewire onto pelvic cavity along with across stricture.

- Insert a cystoscope/ureteroscope to the level of UVJ, and then insert the ureteral catheter through the working channel of cystoscope/ureteroscope. Face ureteral catheter tip onto UVJ orifice. Advance the guide wire across the target stricture where the Stent delivery system to be placed over the guide wire.
- Under fluoroscopic guidance, insert a guide wire across the stricture with sufficient length to allow the Stent delivery system to be placed over the guide wire.
- Free the guidewire from ureteral catheter and cystoscope.
- Ensure the proximal valve of Y-connector locked.
- Flush the inner lumen of delivery system with sterile saline.

④ Stent Deployment Procedure


Figure 2

PRECAUTION: Do not twist delivery system or employ a boring motion during the deployment as this may affect positioning and ultimate function of Stent.

- Following placed guide wire, introduce the delivery system until stricture site.
- Under the fluoroscopic and/or endoscopic guidance, adjust location of stent portion in delivery system to cover stricture length.
- Once the delivery system is considered as in correct position, turn-unlock the Locking knob into counterclockwise direction (2~3 rounds are enough).
- To begin Stent deployment, immobilize the hub in one hand and grasp the Y-connector with the other hand. Gently pull back the Y-connector back along the pusher towards the hub.
- Check fluoroscopic view every second while deployment.


Figure 3

CAUTION Do not push the hub forward or pull backward with the Stent partially deployed. The hub must be securely immobilized. Inadvertent movement of the hub may cause misalignment of the Stent and possible damage to the ureter.

⑤ After Stent Deployment

- Examine the Stent fluoroscopically and/or endoscopically to confirm expansion.
- Carefully remove the introducer system, guidewire and endoscope from the patient. If excessive resistance is felt during removal, wait 3~5 minutes to allow further stent expansion (Place the inner sheath back into the outer sheath as the original state prior to removal.)
- Balloon dilatation inside the stent can be performed if on demand.

11. Perform routine post implant procedure

- Assess the size and structure of the Stent lumen. Full expansion of stent may require 1~2 days.
- Surgeon's experience and discretion can determine the appropriate drug regimen for each patient.
- Observe the patient for development of any complications.

12. Instructions for removal of UVENTA™ Ureteral Stent (see Warnings)

Visually examine the Stent for any tumor over-growth/hyperplasia or encrustation. If the Stent lumen is clear, carefully remove using a forceps. Grasp proximal end of the Stent then carefully retrieve the Stent.

Please note: Since major indication of UVENTA™ Ureteral Stent is for malignant patients, removal procedure is **NOT** post-stent indication. If removal is strongly required due to complications and other reasons, make sure if stent can be safely removed before removal process.

Reuse Precaution Statement

Contents supplied STERILE (ethylene oxide (EO)). Do not use if sterile barrier is damaged. In the event of damaged packaging, call your Taewoong Medical Co., Ltd. representative. UVENTA™ Ureteral Stent is for single use only and only one patient. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and /or cause patient infection or cross infection, including, but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Storage: Store at room temperature (10~40°C).

Disposal Requirements: The introducer system of Uventa™ Ureteral Stent must be properly sealed and disposed in compliance with the regulation of local or hospital at the end of its use.

***Warranty**

Taewoong Medical Co., Ltd. warrants that reasonable care has been used in the design and manufacture of this instrument.

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