


Niti-S Biliary Covered Stent [Flare-Type]

for benign and/or malignant biliary strictures



Features

- The full silicone-covered body prevents tissue-ingrowth
- Flares at both ends prevent migration
- A proximal retrieval string helps for easy removal or repositioning

Ordering Information

Percutaneous				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
TS0604FW	6	4	8.5	50
TS0605FW		5		
TS0606FW		6		
TS0607FW		7		
TS0608FW		8		
TS0609FW		9		
TS0610FW		10		
TS0612FW		12		
TS0804FW		8		
TS0805FW	5			
TS0806FW	6			
TS0807FW	7			
TS0808FW	8			
TS0809FW	9			
TS0810FW	10			
TS0812FW	12			
TS1004FW	10		4	8.5
TS1005FW		5		
TS1006FW		6		
TS1007FW		7		
TS1008FW		8		
TS1009FW		9		
TS1010FW		10		
TS1012FW		12		

Percutaneous, Song Suture				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
TS0804FWS	8	4	9	50
TS0805FWS		5		
TS0806FWS		6		
TS0807FWS		7		
TS0808FWS		8		
TS0809FWS		9		
TS0810FWS		10		
TS0812FWS		12		
TS1004FWS		10		
TS1005FWS	5			
TS1006FWS	6			
TS1007FWS	7			
TS1008FWS	8			
TS1009FWS	9			
TS1010FWS	10			
TS1012FWS	12			

Ordering Information

Endoscopic				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BS0604FW	6	4	8.5	180
BS0605FW		5		
BS0606FW		6		
BS0607FW		7		
BS0608FW		8		
BS0609FW		9		
BS0610FW		10		
BS0612FW		12		
BS0804FW		8		
BS0805FW	5			
BS0806FW	6			
BS0807FW	7			
BS0808FW	8			
BS0809FW	9			
BS0810FW	10			
BS0812FW	12			
BS1004FW	10		4	8.5
BS1005FW		5		
BS1006FW		6		
BS1007FW		7		
BS1008FW		8		
BS1009FW		9		
BS1010FW		10		
BS1012FW		12		

Endoscopic, Short-wire				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BSM0604FW	6	4	8.5	180
BSM0605FW		5		
BSM0606FW		6		
BSM0607FW		7		
BSM0608FW		8		
BSM0609FW		9		
BSM0610FW		10		
BSM0612FW		12		
BSM0804FW		8		
BSM0805FW	5			
BSM0806FW	6			
BSM0807FW	7			
BSM0808FW	8			
BSM0809FW	9			
BSM0810FW	10			
BSM0812FW	12			
BSM1004FW	10		4	8.5
BSM1005FW		5		
BSM1006FW		6		
BSM1007FW		7		
BSM1008FW		8		
BSM1009FW		9		
BSM1010FW		10		
BSM1012FW		12		