


ComVi Biliary Stent

for malignant biliary strictures



Features

- **Triple layered construction:** Biocompatible PTFE membrane is held between inner and outer mesh
 - Unfixed cell structure enables stent to conform to the shape of bile duct
 - PTFE membrane prevents risk of tissue invasion
 - Outer wire mesh prevents the risk of migration
- **Minimum foreshortening** for accurate stent placement
- **Radiopaque marker:** 4 (four) at both covered part ends

ComVi Biliary Stent [Full Covered-Type] for malignant biliary strictures



Ordering Information

Percutaneous				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
TC0604F	6	4	8	50
TC0605F		5		
TC0606F		6		
TC0607F		7		
TC0608F		8		
TC0609F		9		
TC0610F		10		
TC0612F		12		
TC0804F	8	4	8	50
TC0805F		5		
TC0806F		6		
TC0807F		7		
TC0808F		8		
TC0809F		9		
TC0810F		10		
TC0812F		12		
TC1004F	10	4	8	50
TC1005F		5		
TC1006F		6		
TC1007F		7		
TC1008F		8		
TC1009F		9		
TC1010F		10		
TC1012F		12		

ComVi Biliary Stent [Full Covered-Type]

for malignant biliary strictures

Ordering Information

Endoscopic				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BC0604F	6	4	8	180
BC0605F		5		
BC0606F		6		
BC0607F		7		
BC0608F		8		
BC0609F		9		
BC0610F		10		
BC0612F		12		
BC0804F		8		
BC0805F	5			
BC0806F	6			
BC0807F	7			
BC0808F	8			
BC0809F	9			
BC0810F	10			
BC0812F	12			
BC1004F	10		4	8
BC1005F		5		
BC1006F		6		
BC1007F		7		
BC1008F		8		
BC1009F		9		
BC1010F		10		
BC1012F		12		

Endoscopic, Short-wire				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BCM0604F	6	4	8.5	180
BCM0605F		5		
BCM0606F		6		
BCM0607F		7		
BCM0608F		8		
BCM0609F		9		
BCM0610F		10		
BCM0612F		12		
BCM0804F		8		
BCM0805F	5			
BCM0806F	6			
BCM0807F	7			
BCM0808F	8			
BCM0809F	9			
BCM0810F	10			
BCM0812F	12			
BCM1004F	10		4	8.5
BCM1005F		5		
BCM1006F		6		
BCM1007F		7		
BCM1008F		8		
BCM1009F		9		
BCM1010F		10		
BCM1012F		12		

ComVi Biliary Stent [Both Bare-Type] for malignant biliary strictures



Ordering Information

Percutaneous				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
TC0604B	6	4	8	50
TC0605B		5		
TC0606B		6		
TC0607B		7		
TC0608B		8		
TC0609B		9		
TC0610B		10		
TC0612B		12		
TC0804B	8	4	8	50
TC0805B		5		
TC0806B		6		
TC0807B		7		
TC0808B		8		
TC0809B		9		
TC0810B		10		
TC0812B		12		
TC1004B	10	4	8	50
TC1005B		5		
TC1006B		6		
TC1007B		7		
TC1008B		8		
TC1009B		9		
TC1010B		10		
TC1012B		12		

ComVi Biliary Stent [Both Bare-Type]

for malignant biliary strictures

Ordering Information

Endoscopic				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BC0604B	6	4	8	180
BC0605B		5		
BC0606B		6		
BC0607B		7		
BC0608B		8		
BC0609B		9		
BC0610B		10		
BC0612B		12		
BC0804B		8		
BC0805B	5			
BC0806B	6			
BC0807B	7			
BC0808B	8			
BC0809B	9			
BC0810B	10			
BC0812B	12			
BC1004B	10	4	8	180
BC1005B		5		
BC1006B		6		
BC1007B		7		
BC1008B		8		
BC1009B		9		
BC1010B		10		
BC1012B		12		

Endoscopic, Short-wire				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BCM0604B	6	4	8.5	180
BCM0605B		5		
BCM0606B		6		
BCM0607B		7		
BCM0608B		8		
BCM0609B		9		
BCM0610B		10		
BCM0612B		12		
BCM0804B		8		
BCM0805B	5			
BCM0806B	6			
BCM0807B	7			
BCM0808B	8			
BCM0809B	9			
BCM0810B	10			
BCM0812B	12			
BCM1004B	10	4	8.5	180
BCM1005B		5		
BCM1006B		6		
BCM1007B		7		
BCM1008B		8		
BCM1009B		9		
BCM1010B		10		
BCM1012B		12		

ComVi Biliary Stent [End Bare-Type]

for malignant biliary strictures



Ordering Information

Percutaneous				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
TC0804H	8	4	8	50
TC0805H		5		
TC0806H		6		
TC0807H		7		
TC0808H		8		
TC0809H		9		
TC0810H		10		
TC0812H		12		
TC1004H	10	4	8	50
TC1005H		5		
TC1006H		6		
TC1007H		7		
TC1008H		8		
TC1009H		9		
TC1010H		10		
TC1012H		12		

Endoscopic				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BC0804H	8	4	8	180
BC0805H		5		
BC0806H		6		
BC0807H		7		
BC0808H		8		
BC0809H		9		
BC0810H		10		
BC0812H		12		
BC1004H	10	4	8	180
BC1005H		5		
BC1006H		6		
BC1007H		7		
BC1008H		8		
BC1009H		9		
BC1010H		10		
BC1012H		12		