


# ComVi Pyloric/Duodenal Stent [Both Bare-Type]

for intrinsic and/or extrinsic malignant pyloric/duodenal stricture



## Features

- **Triple layered construction:** Biocompatible PTFE membrane tube is held between inner and outer mesh
  - Unfixed cell structure enables stent to conform to tortuous anatomy
  - PTFE membrane prevents the risk of tissue invasion
  - Outer wire mesh prevents the risk of migration
- **Radiopaque marker:** 4 (four) at both covered part ends

# Ordering Information

OTW (Over The Wire)				
Both Bare (5mm)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PC1806B	18	6	14	135
PC1808B		8		
PC1810B		10		
PC1812B		12		
PC2006B	20	6	14	135
PC2008B		8		
PC2010B		10		
PC2012B		12		
PC2206B	22	6	16	135
PC2208B		8		
PC2210B		10		
PC2212B		12		
PC2406B	24	6	16	135
PC2408B		8		
PC2410B		10		
PC2412B		12		
PC2606B	26	6	18	135
PC2608B		8		
PC2610B		10		
PC2612B		12		
PC2806B	28	6	18	135
PC2808B		8		
PC2810B		10		
PC2812B		12		

TTS (Through The Scope)				
Both Bare (5mm)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PCT1806B	18	6	10.5	180
PCT1808B		8		
PCT1810B		10		
PCT1812B		12		
PCT2006B	20	6	10.5	180
PCT2008B		8		
PCT2010B		10		
PCT2012B		12		
PCT2206B	22	6	10.5	180
PCT2208B		8		
PCT2210B		10		
PCT2212B		12		

TTS (Through The Scope)				
Both Bare (5mm)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PCT1806B-22	18	6	10.5	220
PCT1808B-22		8		
PCT1810B-22		10		
PCT1812B-22		12		
PCT2006B-22	20	6	10.5	220
PCT2008B-22		8		
PCT2010B-22		10		
PCT2012B-22		12		
PCT2206B-22	22	6	10.5	220
PCT2208B-22		8		
PCT2210B-22		10		
PCT2212B-22		12		

# Ordering Information

OTW (Over The Wire)				
Both Bare (15mm)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PC1806BA	18	6	14	135
PC1808BA		8		
PC1810BA		10		
PC1812BA		12		
PC2006BA	20	6	14	135
PC2008BA		8		
PC2010BA		10		
PC2012BA		12		
PC2206BA	22	6	16	135
PC2208BA		8		
PC2210BA		10		
PC2212BA		12		
PC2406BA	24	6	16	135
PC2408BA		8		
PC2410BA		10		
PC2412BA		12		
PC2606BA	26	6	18	135
PC2608BA		8		
PC2610BA		10		
PC2612BA		12		
PC2806BA	28	6	18	135
PC2808BA		8		
PC2810BA		10		
PC2812BA		12		

TTS (Through The Scope)				
Both Bare (15mm)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PCT1806BA	18	6	10.5	180
PCT1808BA		8		
PCT1810BA		10		
PCT1812BA		12		
PCT2006BA	20	6	10.5	180
PCT2008BA		8		
PCT2010BA		10		
PCT2012BA		12		
PCT2206BA	22	6	10.5	180
PCT2208BA		8		
PCT2210BA		10		
PCT2212BA		12		

TTS (Through The Scope)				
Both Bare (15mm)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PCT1806BA-22	18	6	10.5	220
PCT1808BA-22		8		
PCT1810BA-22		10		
PCT1812BA-22		12		
PCT2006BA-22	20	6	10.5	220
PCT2008BA-22		8		
PCT2010BA-22		10		
PCT2012BA-22		12		
PCT2206BA-22	22	6	10.5	220
PCT2208BA-22		8		
PCT2210BA-22		10		
PCT2212BA-22		12		