


# ComVi Enteral Colonic Stent [Both Bare-Type]

for intrinsic and/or extrinsic malignant colonic stricture



## Features

- **Triple layered construction: Biocompatible PTFE membrane tube is held between inner and outer mesh**
  - Unfixed cell structure enables stent to conform to tortuous anatomy
  - PTFE membrane prevents the risk of the tissue invasion
  - Outer wire mesh prevents the risk of migration
- **Radiopaque marker: 4 (four) at both covered part ends**

## Ordering Information

TTS (Through The Scope)				
5mm bare				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CCT1806B	18	6	10.5	220
CCT1808B		8		
CCT1810B		10		
CCT1812B		12		
CCT2006B	20	6	10.5	220
CCT2008B		8		
CCT2010B		10		
CCT2012B		12		
CCT2206B	22	6	10.5	220
CCT2208B		8		
CCT2210B		10		
CCT2212B		12		

TTS (Through The Scope)				
15mm bare				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CCT1806BA	18	6	10.5	220
CCT1808BA		8		
CCT1810BA		10		
CCT1812BA		12		
CCT2006BA	20	6	10.5	220
CCT2008BA		8		
CCT2010BA		10		
CCT2012BA		12		
CCT2206BA	22	6	10.5	220
CCT2208BA		8		
CCT2210BA		10		
CCT2212BA		12		

# Ordering Information

OTW (Over The Wire)				
5mm bare				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CC1806B	18	6	14	70
CC1808B		8		
CC1810B		10		
CC1812B		12		
CC2006B	20	6	14	70
CC2008B		8		
CC2010B		10		
CC2012B		12		
CC2206B	22	6	16	70
CC2208B		8		
CC2210B		10		
CC2212B		12		
CC2406B	24	6	16	70
CC2408B		8		
CC2410B		10		
CC2412B		12		
CC2606B	26	6	18	70
CC2608B		8		
CC2610B		10		
CC2612B		12		
CC2806B	28	6	18	70
CC2808B		8		
CC2810B		10		
CC2812B		12		
CC3006B	30	6	18	70
CC3008B		8		
CC3010B		10		
CC3012B		12		

OTW (Over The Wire)				
15mm bare				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CC1806BA	18	6	14	70
CC1808BA		8		
CC1810BA		10		
CC1812BA		12		
CC2006BA	20	6	14	70
CC2008BA		8		
CC2010BA		10		
CC2012BA		12		
CC2206BA	22	6	16	70
CC2208BA		8		
CC2210BA		10		
CC2212BA		12		
CC2406BA	24	6	16	70
CC2408BA		8		
CC2410BA		10		
CC2412BA		12		
CC2606BA	26	6	18	70
CC2608BA		8		
CC2610BA		10		
CC2612BA		12		
CC2806BA	28	6	18	70
CC2808BA		8		
CC2810BA		10		
CC2812BA		12		
CC3006BA	30	6	18	70
CC3008BA		8		
CC3010BA		10		
CC3012BA		12		