


Niti-S Enteral Colonic Covered Stent



Features

- **Fixed cell with braided construction**
- **Silicone covering and soft round ends**
- **Visible green suture for easy removal**
- **Radiopaque marker: 4 (four) at both ends & 2 (two) in the middle**

Niti-S Enteral Colonic Covered Stent [Fully Covered-Type]

for intrinsic and/or extrinsic malignant and/or benign colonic stricture

Ordering Information

TTS (Through The Scope)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CST1806F	18	6	10.5	220
CST1808F		8		
CST1810F		10		
CST1812F		12		
CST1814F		14		
CST1815F		15		
CST2006F	20	6	10.5	220
CST2008F		8		
CST2010F		10		
CST2012F		12		
CST2014F		14		
CST2015F		15		
CST2206F	22	6	10.5	220
CST2208F		8		
CST2210F		10		
CST2212F		12		
CST2214F		14		
CST2215F		15		

OTW (Over The Wire)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CS2006F	20	6	20	70
CS2008F		8		
CS2010F		10		
CS2012F		12		
CS2014F		14		
CS2015F		15		
CS2016F	22	16	20	70
CS2206F		6		
CS2208F		8		
CS2210F		10		
CS2212F		12		
CS2214F		14		
CS2215F	24	15	22	70
CS2216F		16		
CS2218F		18		
CS2223F		23		
CS2406F		6		
CS2408F		8		
CS2410F	24	10	22	70
CS2412F		12		
CS2414F		14		
CS2415F		15		
CS2416F		16		
CS2418F		18		
CS2423F	26	23	22	70
CS2606F		6		
CS2608F		8		
CS2610F		10		
CS2612F		12		
CS2614F		14		
CS2615F	28	15	22	70
CS2616F		16		
CS2806F		6		
CS2808F		8		
CS2810F		10		
CS2812F		12		
CS2814F	28	14	22	70
CS2815F		15		
CS2816F		16		

OTW (Over The Wire)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CS1806F	18	6	16	70
CS1808F		8		
CS1810F		10		
CS1812F		12		
CS1814F		14		
CS1815F		15		
CS1816F		16		

Niti-S Enteral Colonic Covered Stent [Both Bare-Type]

for intrinsic and/or extrinsic malignant colonic stricture

Ordering Information

TTS (Through The Scope)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CST1806B	18	6	10.5	220
CST1808B		8		
CST1810B		10		
CST1812B		12		
CST1814B		14		
CST1815B		15		
CST2006B	20	6	10.5	220
CST2008B		8		
CST2010B		10		
CST2012B		12		
CST2014B		14		
CST2015B		15		
CST2206B	22	6	10.5	220
CST2208B		8		
CST2210B		10		
CST2212B		12		
CST2214B		14		
CST2215B		15		

OTW (Over The Wire)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CS2006B	20	6	20	70
CS2008B		8		
CS2010B		10		
CS2012B		12		
CS2014B		14		
CS2015B		15		
CS2016B	22	16	20	70
CS2206B		6		
CS2208B		8		
CS2210B		10		
CS2212B		12		
CS2214B		14		
CS2215B	24	15	22	70
CS2216B		16		
CS2406B		6		
CS2408B		8		
CS2410B		10		
CS2412B		12		
CS2414B	26	14	22	70
CS2415B		15		
CS2416B		16		
CS2606B		6		
CS2608B		8		
CS2610B		10		
CS2612B	28	12	22	70
CS2614B		14		
CS2615B		15		
CS2616B		16		
CS2806B		6		
CS2808B		8		
CS2810B	10			
CS2812B	28	12	22	70
CS2814B		14		
CS2815B		15		
CS2816B		16		

OTW (Over The Wire)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CS1806B	18	6	16	70
CS1808B		8		
CS1810B		10		
CS1812B		12		
CS1814B		14		
CS1815B		15		
CS1816B		16		

Niti-S Enteral Colonic Covered Stent [End Bare-Type]

for intrinsic and/or extrinsic malignant colonic stricture

Ordering Information

TTS (Through The Scope)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CST1806H	18	6	10.5	220
CST1808H		8		
CST1810H		10		
CST1812H		12		
CST1814H		14		
CST1815H		15		
CST2006H	20	6	10.5	220
CST2008H		8		
CST2010H		10		
CST2012H		12		
CST2014H		14		
CST2015H		15		
CST2206H	22	6	10.5	220
CST2208H		8		
CST2210H		10		
CST2212H		12		
CST2214H		14		
CST2215H		15		

OTW (Over The Wire)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CS2006H	20	6	20	70
CS2008H		8		
CS2010H		10		
CS2012H		12		
CS2014H		14		
CS2015H		15		
CS2016H	22	6	20	70
CS2208H		8		
CS2210H		10		
CS2212H		12		
CS2214H		14		
CS2215H		15		
CS2216H	24	6	22	70
CS2408H		8		
CS2410H		10		
CS2412H		12		
CS2414H		14		
CS2415H		15		
CS2416H	26	6	22	70
CS2608H		8		
CS2610H		10		
CS2612H		12		
CS2614H		14		
CS2615H		15		
CS2616H	28	6	22	70
CS2808H		8		
CS2810H		10		
CS2812H		12		
CS2814H		14		
CS2815H		15		
CS2816H	18	6	16	70
CS1808H		8		
CS1810H		10		
CS1812H		12		
CS1814H		14		
CS1815H		15		
CS1816H	16			



OTW (Over The Wire)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CS1806H	18	6	16	70
CS1808H		8		
CS1810H		10		
CS1812H		12		
CS1814H		14		
CS1815H		15		
CS1816H	16			